

# BPI Super Cleaner™

For use only by qualified personnel in a laboratory environment.

Due to high operating temperature, access should be restricted.

Also available for the BPI Super Ultrasonic Cleaner™

IMMERSION TRAY	
BPI# 67300	Immersion tray
LID	
BPI# 67500	Stainless steel lid
BEAKER RACK	
BPI# 67100	Holds one 600 ml beaker
BEAKERS & ACCESSORIES	
BPI# 62300	600 ml beaker
BPI# 63000	Lid for 600 ml beaker
CONTACT LENS RACKS & VIALS	
BPI# 67800	Contact lens rack with 12 vials
BPI# 362100	Vial for contact lens rack

## CAUTION

DO NOT IMMERSE UNIT. DO NOT ALLOW ANY LIQUID INTO INTERIOR OF UNIT. This can void your warranty.



BPI Super Cleaner™

## About Your BPI Cleaner

The BPI Super Cleaner (BPI #7500) is a dual axis ultrasonic that has cleaning power pulses which provide random-sized cavitation bubbles with no standing waves. It has two ON-OFF switches (one for each axis). The Horizontal Axis projects ultrasound waves from back to front of the system; the Vertical Axis projects from the bottom to the top. For general cleaning, both axes may be used for speed and efficiency.

The Cleaner will give best results for general cleaning when used with BPI Ultrasonic "C", which has special chemical abilities that sequester ionic residue and may be used in most all cleaning applications. DO NOT USE ULTRASONIC "C" FOR CONTACT LENS CLEANING.

For optimum cleaning, it is most important to fill the tank to the correct level to ensure proper frequency resonance. An optional immersion tray is available for general parts cleaning.

## Specifications

Your Cleaner has a non-removable stainless steel tank contained in a stainless steel chassis. The Cleaner requires 115 volts, 60 Hz unless otherwise ordered. It is also available by special order in 220 volts.

HEIGHT	WIDTH	LENGTH	VOLTAGE	WEIGHT	FUSE	AMPERAGE
7.5 in.	11.25 in.	10.5 in.	115 v. or 220 v.	11 lbs	2 amp (Two 2 amp fuses)	2 amp
191 mm	285 mm	266 mm		5 kg		
TANK DEPTH	TANK WIDTH	TANK LENGTH	TANK CAPACITY	PEAK POWER	SET-UP KIT	
4 in.	9.5 in.	5.5 in.	3 Quarts	250 watts	1 Lid	
102 mm	241 mm	140 mm	2839 ml		1 Instruction manual	

## Unpacking

When unpacking your ultrasonic cleaner, please check to ensure that no concealed damage occurred in transit. If such is noted, save the shipping carton and immediately notify the shipping company's damage control inspector in your area so a claim may be processed. Failure to do this may void any future claim and replacement. Also, call BPI Customer Service so arrangements for a replacement may be made.

## Notes

- NEVER turn the ultrasonic on when tank is empty as PERMANENT DAMAGE may result. The tank should be at least ¾ full of water or diluted Ultrasonic "C" solution when the instrument is on.
- Always UNPLUG unit from power source before filling tank.
- DO NOT immerse hands in cleaning solution while unit is in operation.
- NEVER remove ground prong from power plug as it ensures safe electrical operation.
- Ultrasonic energy generates heat in the solution. You should periodically check the item(s) being cleaned to avoid unnecessary thermal exposure. Use caution when immersing hand.

## General Cleaning

1. Fill tank with tap water approximately ¾" from the top; wipe excess water from the cabinet and work area and plug in the power cord.
2. Add 3 capfuls of BPI Ultrasonic "C" cleaning concentrate to water in tank.
3. Place items to be cleaned in tank. When cleaning frames with lenses, make sure to face the lenses upward so they do not rest on the bottom of the tank.
4. Turn on both the Horizontal and Vertical Axes.
5. For most cleaning operations, 1-3 minutes is sufficient; more difficult jobs will require

proportionately longer times. You will notice that the Cleaner will not harm the thinnest gold platings.

6. When the cleaning is completed, turn off both power switches; remove cleaned items promptly.

7. If cleaning jewelry, DO NOT place pearls or opals in the Cleaner. While other gems may be cleaned in the Cleaner, BPI assumes no liability for these items.

## Soft Contact Lenses

The cleaning of a soft contact lens is a difficult process even with the best of chemicals. The BPI Cleaner will remove most mucoprotein deposits by fracturing them into microscopic particles, breaking their structure so they can no longer attach themselves to the lens matrix. The BPI Cleaner changes the frequency of oscillation and phase angle of the wave fronts in a continuous manner that produces random-sized bubbles that fracture the protein and mucous matrix. The yellow color in the soft lens due to thermal modification of the material is a physical part of the lens and cannot be removed.

## Cleaning Soft Contact Lenses

1. Place each soft lens in the special bottle provided with the C.L. Rack; the bottle caps are marked R and L (right and left). DO NOT use the nylon mesh insert with soft contact lenses.
2. Place the lens manufacturer's suggested cleaning or storage solution in bottle, covering the lens by ½" of fluid. Screw lid on tightly; place bottle in special rack.
3. Turn both Horizontal and Vertical switches on for 15-20 minutes.
4. When lenses have finished cleaning cycle, turn bottle upside-down and shake solution and lens to the cap; loosen cap and let excess fluid run out around bottle and lid so lens may be easily removed.
5. Follow lens manufacturer's suggested sterilization process before placement into patient's eyes.

## Cleaning Hard Contact Lenses

1. First, be sure to place the nylon mesh insert in the bottom of the bottle before putting the hard contact lens in it. This insert is furnished to prevent the formation of a vacuum which may cause the contact lens to adhere to the bottom of the bottle. Removal of a lens may then be extremely difficult and cause damage to the lens.
2. Insert lens and add just enough water to cover it; then add 15 drops of C.L. Polish to the bottle.
3. Turn both Horizontal and Vertical Axes on for 15-20 minutes.
4. Remove lens from bottle; wash thoroughly with water to remove polishing compound.
5. Then clean with contact lens cleaner recommended by the lens manufacturer to remove any C.L. Polish that may be remaining.

## C.L. Polish

BPI C.L. Polish (for HARD CONTACT LENSES ONLY) is a combination of selected polishing particles which have uniform shape, size, and hardness. These particles have been suspended in a manner to give maximum heat dispersion and transfer while providing a maximum number of polishing particles at the lens surface at any given time.

**Precautions to be observed when using C.L. Polish:** IF FLUID COMES INTO CONTACT WITH EYES, WASH IMMEDIATELY WITH WATER. IF IRRITATION PERSISTS, CONTACT PHYSICIAN.

## Questions? Ordering...

If you have any questions about the use of your Ultrasonic Cleaner or need technical product information, call the BPI Research Department. For information about any other BPI product and to order supplies, please give us a toll-free call at the number for your area.

© 2000 BPI. All specific product names mentioned herein are trademarks of Brain Power Incorporated, Miami, Florida, USA. (Unless otherwise stated), BPI is a registered trademark with the US Patent Office and with similar offices in other countries. MANUAL FILE# M2065

BPI Super Cleaner™

BPI# 7501 (115V)  
BPI# 207501 (220V)